

- Motorized recline
- Hydraulic low lift 9" (22.9cm) stroke with maximum seat height of 29" (73.7cm) or hydrualic high lift 12.5" (31.8cm) stroke with maximum seat height of 35" (88.9cm)
- User programmable memory two position control
- Illuminated membrane switch control on both chair sides
- Solid cast aluminum frame, armrest, footrest, and base
- One hand articulating headrest
- Securely locking hand brake with handle on both chair sides
- Footrest



Technical Specifications

Elevation Stroke	Low: 9" (22.9cm) / High: 12.5" (31.8cm)
Voltage & Frequency	115V; 5A; 60Hz
	220V; 3A; 50/60Hz
Rotation Angle	355°
Chair Capacity	500lb (226.8kg)
Shipping Weight	286lb (129.7kg)
ETL Listed	Complies with UL and CSA UL-60601-1 CAN/CSA - C22-2 No. 60-1-M90

Upholstery Selections


Standard Colors



Designer Series Colors







Topcon Medical Systems, Inc.

111 Bauer Drive, Oakland, NJ 07436, U.S.A. t 800.223.1130 f 201.599.5248 • www.topconmedical.com © 2014 Topcon Medical Systems, Inc. MCA# 713