



500 Series Hydraulic Surgical Stools

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



500 Series Hydraulic Surgical Stools

500 Series Surgical Stools

The base exceeds the circumference of the stool to provide stability, even in the sit-stand position. Swivel of stool seat is locked by lifting floor pedal. A convenient foot pedal provides hands-free adjustment of stool height.

Choice of Three Heights

Standard: 17½" to 25", medium: 20½ to 28", high: 24" to 31½".

Finished for Durability

- Select from 14, easy-to-maintain vinyl colors. Other vinyl colors available by request.
- All exposed metal, with exception of the aluminum base, is chrome-plated steel.
- Base and legs finished in tough, black epoxy enamel.

Optional Foot Ring

Order a foot ring with your stool, or at a later date.

Exceptional Stability

All 500 series stools have five legs to evenly distribute weight and ensure stability.

Tech Specs

Base diameter 19¾"; footrest: 16" diameter; 2" spacer available to elevate ring 2" to 4".



	540	546	548	556
Seat adjustment	Foot Pedal	Foot Pedal	Foot Pedal	Foot Pedal
Seat	15" circular	15" circular	14 x 16" rectangular	15" circular
Backrest – tilt, horizontal and vertical adjustment	None	Yes	Yes	Yes
Wrap-around arm	No	No	No	Yes
Seat	2" polyfoam	2" polyfoam	2" polyfoam	2" polyfoam
Five legs for extra stability	Yes	Yes	Yes	Yes
Foot pedal	Yes	Yes	Yes	Yes
Floor lock	Optional	Optional	Optional	Optional
Shipping weight	47 lbs	51 lbs	57 lbs	63 lbs