



BA 904 PORTABLE SLIT LAMP

High-quality examinations, anywhere...

02 | 03 POWERFUL, HAND-HELD SLIT LAMP

BA 904

Haag-Streit optics included

The BA 904 is the only hand-held slit lamp that boasts Haag-Streit's unsurpassed optical system.

The robust, ergonomically-designed portable BA 904 slit lamp can provide high-quality examinations, even when you practice in more than one consultancy room, or perform field-based eye examinations.

The BA 904 can be used for hand-held operation to examine patients who cannot comfortably sit at a larger slit lamp, yet can be easily reassembled for more traditional joystick/headrest operation. Whichever method is used, the portability does not diminish the quality of the examination.

Both reliable and accurate, the BA 904 offers a Dioptic range of -8 to +8, an interpupillary range of 53 to 95. It offers a wider field of view, compared to other portable slit lamps available.

Both light and portable, the complete kit fits into the specially-designed light-weight suitcase. The integrated rechargeable NICAD battery allows up to 45 minutes of continuous use per battery (two batteries supplied).

High-quality Haag-Streit optics

The BA 904 is the only portable slit lamp that incorporates World-class Haag-Streit optics. These superior optics offer high-resolution examinations; boasting 10x objectives, as standard and 16x objectives, as an option.

Versatile

A dual-function slit lamp, the BA 904's ergonomic design enables both hand-held examinations and traditional-style examinations, using the easy-to-assemble head and chin rest.

Portable

Extremely light and portable, the BA 904 effortlessly fits into a specially-designed suitcase, which is supplied as standard. It is easy to transport between consultancy rooms or field-based locations.

Safe & comfortable

A well-balanced instrument, the BA 904 portable, hand-held slit lamp is both safe and comfortable for the patient. It filters out hazardous radiation, which means less damage to the retina when using a non-contact lens.

Integrated battery

Its integrated, rechargeable NiCad battery allows up to 90 minutes examination time, away from the consultancy room. Two batteries are supplied, as standard.

Two BA 904 models available

You can choose the portable slit lamp to suit your individual requirements from the two models that are available;

Set A

BA 904 portable slit lamp, including; head and chin rest stand, two energy packs, charger, power supply and large case (Part Number: 5606011).

Set B

BA 904 portable slit lamp, including; two energy packs, charger, power supply, parking unit and small case (Part Number: 5606013).

Meets IEC 60601-1 3rd Edition & regulatory standards

The BA 904 has been updated to meet the stringent IEC 60601-1 3rd edition and regulatory standards. It has undergone rigorous safety, Electromagnetic Compatibility (EMC) and environmental testing, including; temperature range, vibration and humidity testing, as-well-as minor design improvements. All plastic parts are completely flame-retardant.



HAAG-STREIT UK

Edinburgh Way
Harlow
Essex
CM20 2TT
Phone (01279) 456261
info@haag-streit-uk.com
www.haagstreituk.com

Members of HAAG-STREIT Group

HAAG-STREIT Holding AG
www.haag-streit-holding.com

HAAG-STREIT Deutschland GmbH
www.haag-streit.de

HAAG-STREIT AG
www.haag-streit.com

IPro GmbH
www.ipro.com

SPECTROS AG
www.spectros.ch

CLEMENT CLARKE Ltd.
www.clement-clarke.com

HAAG-STREIT Medtech AG
www.haag-streit-medtech.com

HAAG-STREIT UK
www.haag-streit-uk.com

HAAG-STREIT France EURL
www.haag-streit.fr

John Weiss Ltd.
www.johnweiss.com

HAAG-STREIT Far East
www.haag-streit-fareast.com

HAAG-STREIT USA
www.haag-streit-usa.com

HAAG-STREIT Surgical GmbH
www.haag-streit-surgical.com

Reliance Medical Inc.
www.haag-streit-usa.com

Möller-Wedel GmbH & Co KG
www.moeller-wedel.com

Asetronics AG
www.asetronics.ch

Möller-Wedel Optical GmbH
www.moeller-wedel-optical.com

ComLab AG
www.comlab.ch